

Ormco[™]
Your Practice. Our Priority.

The Mini Diamond System



Mini Diamond®

The Mini Diamond reflects years of research and development. The goal was to improve aesthetics by miniaturizing the bracket with no compromise in function.

Mini Diamonds offer the same benefits in regard to accurate bracket placement as their traditional counterparts plus:


- Improved aesthetics and hygiene with a 30% reduction in overall bracket size.
- Improved patient comfort with enhanced facial contours.
- Improved utility with refined Hemi Hooks™ and permanent I.D. marking system.
- Confidence in performance. The Mini Diamond Bracket is made from 17-4 stainless steel, the hardest steel brackets can be manufactured from
- Available in a variety of prescriptions including MBT*
- Gingival offset pads if required



Mini Diamond®

Mini Diamond Marking System

- Open dimple (°) = Maxillary
- Red dimple (•) = Mandibular
- On Mini Diamond or other angulated assemblies, the dimple is on the distogingival tie-wing (twin) or gingival tie-wing (single).
- On non-angulated assemblies, dimples can be on either gingival tie-wing.
- On lower bicuspid, where torques differ, the 2nd bicuspid have dimples on both gingival tie-wings.
- The lower cuspids and bicuspid also feature a scribe line on the occlusal tie-wings.



Item	Description	Item	Description
351-0130	Upper Central Right +I4T+5A 22 10/BOX	351-0025	Lower Anteriors -1T 0A 022 10/BOX
351-0131	Upper Central Left +I4T+5A 22 10/BOX	351-0035	Lower Anteriors -6T 0A 022 10/BOX
351-0140	Upper Central Right +22T+5A 22 10/BOX	351-1422	Lower Cuspid Right W/HK 0T +6A .022 10/BOX
351-0141	Upper Central Left +22T+5A 22 10/BOX	351-1423	Lower Cuspid Left W/HK 0T +6A .022 10/BOX
351-0250	Upper Lateral Right +7T +8A 022 10/BOX	351-1432	Lower Cuspid Right W/HK -7T +6A 022 10/BOX
351-0251	Upper Lateral Left +7T +8A 022 10/BOX	351-1433	Lower Cuspid Left W/HK -7T +6A 022 10/BOX
351-1320	Upper Cuspid Right W/HK -7T +10A .022 10/BOX	351-0615	Lower 1st Bicuspid Universal -11T 0A 022 10/BOX
351-1321	Upper Cuspid Left W/HK -7T +10A .022 10/BOX	351-0635	Lower 1st Bicuspid Universal G/O -11T 0A 022 10/BOX
351-1330	Upper Cuspid Right W/HK 0T+10A 022 10/BOX	351-0725	Lower 2nd Bicuspid Universal -17T 0A 022 10/BOX
351-1331	Upper Cuspid Left W/HK 0T+10A 022 10/BOX	351-0745	Lower 2nd Bicuspid Universal G/O -17T 0A 022 10/BOX
351-0514	Upper Bicuspid Universal -7T 0A 022 10/BOX		
351-0534	Upper Bicuspid Universal G/O -7T 0A 022 10/BOX		

*available with non hooked cuspids and hooks on all premolars
* not endorsed by Drs. McLaughlin, Bennett, or Trevisi.*

ACCENT^{mini}



- Contoured, low profile design for a precision fit with less occlusal interference and greater patient comfort
- Single-piece with patented laser-etched base for optimal bond strength and reliability
- Scribe line for accurate bracket placement and precision finishing
- Flared mesial slot opening for easy wire changes
- Bendable “J” hook for greater elastic retention and patient comfort
- Scaler notch along the M-D length for easy adhesive loading and angulation adjustments
- Available in MBT* and Damon System prescriptions

Damon® Rx

Upper	Torque	Angulation	Distal Offset	.022 L	.022 R
First Molar	-18	0	12	513 0163	513 0162
Second Molar	-27	0	+6	513 0171	513 0170
Lower	Torque	Angulation	Distal Offset	.022 L	.022 R
First Molar	-28	2	2	513 0165	513 0164
Second Molar	-10	0	+5	513 0173	513 0172

MBT Rx*

Upper	Torque	Angulation	Distal Offset	.018 L	.018 R	.022 L	.022 R
First Molar	-14	0	10	340-1365	340-1364	342-1365	342-1364
Second Molar	-14	0	10	340-1375	340-1374	342-1375	342-1374
Lower	Torque	Angulation	Distal Offset	.018 L	.018 R	.022 L	.022 R
First Molar	-20	0	0	340-1367	340-1366	342-1367	342-1366
Second Molar	-10	0	5	340-1377	340-1376	342-1377	342-1376



Upper 2nd Molar Tube



Lower 2nd Molar Tube



Lower 1st Molar Tube




Upper 1st Molar Tube

Archwires

To complete the total conventional offering, Ormco are pleased to offer affordable yet high quality Niti and stainless steel wires for all your clinical needs

Part Number	Description
205-0101	+ NITI 2000 Truarch .014 DIA Upper Med BP/100
205-0102	+ NITI 2000 Truarch .016 DIA Upper Med BP/100
205-0103	+ NITI 2000 Truarch .018 DIA Upper Med BP/100
205-0111	+ NITI 2000 Truarch .014 DIA Lower Med BP/100
205-0112	+ NITI 2000 Truarch .016 DIA Lower Med BP/100
205-0113	+ NITI 2000 Truarch .018 DIA Lower Med BP/100
211-0521	+ NITI 2000 Truarch .016X.022 Lower Med BP/100
211-0523	+ NITI 2000 Truarch .017X.025 Lower Med BP/100
211-0524	+ NITI 2000 Truarch .019X.025 Lower Med BP/100

Part Number	Description
227-9003	+ TRUARCH SS LM 017X025 ND PK/10
227-9031	+ TRUARCH SS UM 017X025 ND PK/10
227-9034	+ TRUARCH SS UM 019X025 ND PK/10
227-9006	+ TRUARCH SS LM 019X025 ND PK/10



Ortho Solo™

- Fluoride-releasing universal sealant
- No air-dry or light-cure step is needed
- 24% filled – giving increased bond strength
- Displaces moisture
- Also available in sticks - consistent and accurate dispensing, minimizing waste & reduces the risk of cross contamination

Ortho Solo Kit 740-0270. 4 x 5ml bottles & brushes
 Ortho Solo Refill 740-0271. 1 x 5ml bottle
 Box of 10 Sticks (1.2 ml total) 740-0269

