

TUBES

Accent Tubes

BONDABLE ACCENT TUBES

The Molar Accent tubes exist in single or double tube for the first molar and only in single tube for the second molar, available in .018" x .025" and .022" x .028". These tubes are non convertible, low profile with no tie wings limiting occlusal interference.

- ▶ Made of 17-4 grade stainless steel.
- ▶ Conical mesial opening for easy archwire insertion.
- ▶ Available with distal extension on second molars.

Single



Double



WELDABLE ACCENT TUBES

- ▶ Made of 17-4 grade stainless steel.
- ▶ Conical mesial opening for easy archwire insertion.
- ▶ Available with distal extension on second molars.

Single



Double

