

TUBES

Accent Mini Tubes

ACCENT MINI TUBES

Latest addition to the Ormco tubes range, Accent mini molar tubes feature a “low profile” reducing occlusal interference, providing greater patient comfort and facilitating positioning during bonding. Manufactured with 17-4 grade stainless steel using metal injection molding (MIM); laser-etched pad for optimal bonding.

- ▶ Adjustable «T» hook for the 1st lower molars to facilitate cross-elastic hold.
- ▶ Funneled mesial slot opening for easy wire insertion.
- ▶ Occlusal-gingival notch for easy adhesive loading and bracket placement.
- ▶ Horizontal marking for easy tube positioning.

1st upper molar



1st lower molar



2nd upper molar



2nd lower molar

