




Available in Roth Type, MBT*, .020" x .028",
Extra Torque.
Single patient kits, or in bags of 10.

Chic.



Made in Japan by



So Chic!

Visit www.gcorthodontics.eu for further information.



Printed on a 60 % recycled paper - *MBT is a trademark of 3M Unitek - 99-14E0-1100 Rev04 - 12/2015 - *Roth is a registered trademark of Roth Licensing LLC.

GC Orthodontics
GLOBAL EXPERTS IN ORTHODONTICS

Chic.

Hi Aesthetic

Chic. is a ceramic bracket made from highly translucent polycrystalline alumina and features a 100% mechanical base for efficient bonding with any type of adhesive. Created from a patented firing process, the bracket provides unrivalled smoothness for increased patient comfort. The long axis scribe line dissolves with brushing.

Initialloy and BIO-ACTIVE Aesthetic Archwires

Combine Chic. brackets with our rhodium coated Initialloy and BIO-ACTIVE archwires for outstanding aesthetics.



Rhomboid shaped bracket with axis line for effortless placement, non-permanent line disappears after brushing


Disto-gingival color ID

Chamfered slot for easy wire insertion

Highly translucent polycrystalline allumine bracket

True twin bracket design for an uncompromised treatment

Chic.



Torque in base



Wide area under tie-wing to facilitate ligation

100% mechanical base (identical to EXPERIENCE C base)

Color Codes



Disto-gingival color code aides with bracket identification in order to prevent bonding errors.

