

BRACKETS

Lingual Self-Ligating

ALIAS BRACKET

Alias, the first passive self-ligating square slot lingual straight-wire bracket developed through collaboration with world-renowned lingual leaders, Drs Kyoto Takemoto and Giuseppe Scuzzo, offers an aesthetic alternative to buccal brackets or aligners with the efficiency of Straight-Wire™ mechanics. Alias' unique square slot allows for optimized movements and greater interbracket distances make the aligning stage easier. Easier bonding, especially on small maloccluded incisors, the technical features of Alias facilitate lingual orthodontics and enhance clinical performance.



ALIAS
LINGUAL STRAIGHTWIRE BRACKET SYSTEM
Developed with Drs. Kyoto Takemoto and Giuseppe Scuzzo

► Innovation

A patented square slot .018" x .018" provides optimal control and results.

► Comfort

Greater patient comfort with slim, rounded brackets.

► Efficacy

Simplified Straight-wire mechanics.

► Simplicity

Easy opening and closing thus decreasing chair time.



Tooth	Torque	Angulation	Distal Offset	Hook	Slot Size	
					.018 L	.018 R
MAXILLARY						
Anteriors	60 °	0 °	0 °		369-0020	369-0020
Cuspid	60 °	8 °	+06 °		369-0031	369-0030
2nd Bicuspid	0 °	0 °	0 °		369-0050	369-0050
1st Molar	10 °	0 °	+05 °		369-0061	369-0060
2nd Molar	10 °	0 °	+05 °	Yes	369-1071	369-1070
MANDIBULAR						
Anteriors	45 °	0 °	0 °		369-0002	369-0002
Cuspid	45 °	8 °	+06 °		369-0011	369-0012
1st Molar	0 °	0 °	0 °		369-0080	369-0080
2nd Molar	0 °	0 °	0 °	Yes	369-1091	369-1090
MAXILLARY/MANDIBULAR						
Bicuspid	0 °	0 °	0 °		369-0040	369-0040
7-7/7-7				Yes	740-0430	740-0430