

ADHESIVES, CEMENTS & PRIMERS

Composites

BLÜGLOO™

Blüglöo™ light-curing color contrast adhesive. The color contrast of Blüglöo™ during bonding helps facilitate fast, accurate bracket placement and makes it easy to see unwanted flash. As Blüglöo™ warms to body temperature, the color virtually disappears, remaining clear throughout treatment. When debonding, simply introduce a short blast of cool air or water to lower the bonding surface temperature, and Blüglöo™ turns blue again for thorough cleanup.

- ▶ Color on demand: blue at bonding, clear at treatment, and blue at cleanup.
- ▶ Made for aesthetics.
- ▶ Easy to use.
- ▶ Generous working time.



Designation	Kit Content	Part Number
Blüglöo™ Syringe Kit	2 syringes (2x4 g) Ortho Solo (5 ml) Spatulas (75) Etching solution (10 ml) Applicator brushes (38)	740-0274
Blüglöo™ Unidose Kit	4 packs of 10 unidoses (2x4g) Ortho Solo (5 ml) Spatulas (75) Etching solution (10 ml) Applicator brushes (38) Unidose dispensing gun	740-0275

Refill	Part Number
Blüglöo syringe (4 g)	740-0272
2 packs of 10 Blüglöo unidoses (4 g in total)	740-0273
1 Ortho Solo bottle (5 ml)	740-0271
Bottle of etching solution (10 ml)	740-0038
Brushes (75)	740-0041
Unidose dispensing gun	740-0200
Spatulas (75)	740-0075