

Comfort Line 3D buccal tubes


Info
N° 216

Comfort Line **3D**

Comfort Line 3D buccal tubes

Comfort Line 3D buccal tubes are the newest generation of our proven classics providing optimised bonding strength via an anatomically shaped, smaller pad. These MIM-manufactured one piece buccal tubes have all advantages of Comfort Line buccal tubes and are also nickel free. The smart design of the 3D-base provides optimal bonding strength whilst maintaining safe and comfortable debonding.

anatomical indent



Comfort Line 3D buccal tubes are also available with an anatomical indent for precision placement.



Maxillary	Tooth	Torque	Angulation	Distal offset	Slot				
						without anatomical indent		with anatomical indent	
						right	left	right	left
	76 67	0°	—	0°	.018" x .028"	748M8211	748M8311	748K8211	748K8311
					.022" x .028"	748M8221	748M8321	748K8221	748K8321
	76 67	-14°	—	12°	.018" x .028"	728M0804	727M0804	728K0804	727K0804
					.022" x .028"	728M0802	727M0802	728K0802	727K0802
	76 67	-10°	—	12°	.018" x .028"	774M0804	773M0804	774K0804	773K0804
					.022" x .028"	774M0802	773M0802	774K0802	773K0802
	76 67	-14°	—	14°	.018" x .028"	741M0804	740M0804	741K0804	740K0804
					.022" x .028"	741M0802	740M0802	741K0802	740K0802

Mandibular

	76 67	0°	—	0°	.018" x .028"	748M8311	748M8211	748K8311	748K8211
					.022" x .028"	748M8321	748M8221	748K8321	748K8221
	76 67	-10°	—	0°	.018" x .028"	728M1804	727M1804	728K1804	727K1804
					.022" x .028"	728M1802	727M1802	728K1802	727K1802
	76 67	-20°	—	0°	.018" x .028"	728M1714	727M1714	728K1714	727K1714
					.022" x .028"	728M1712	727M1712	728K1712	727K1712
	76 67	-25°	—	12°	.018" x .028"	774M1804	773M1804	774K1804	773K1804
					.022" x .028"	774M1802	773M1802	774K1802	773K1802
	76 67	-30°	—	0°	.018" x .028"	748M8315	748M8215	748K8315	748K8215
					.022" x .028"	748M8325	748M8225	748K8325	748K8225

Buccal tubes

Roth Prescription



Maxillary	Tooth	Torque	Angulation	Distal offset	Slot	for welding		on small pad	
						right	left	right	left
 occlusal headgear tube convertible	6 6	-14°	-	14°	.018" x .028"	741-0703	740-0703	-	-
					.022" x .028"	741-0701	740-0701	-	-
 convertible	6 6	-14°	-	14°	.018" x .028"	741-0727	740-0727	-	-
					.022" x .028"	741-0725	740-0725	-	-
 gingival headgear tube convertible	6 6	-14°	-	14°	.018" x .028"	741-0713	740-0713	-	-
					.022" x .028"	741-0711	740-0711	-	-
 occlusal headgear tube convertible	6 6	-14°	-	14°	.018" x .028"	741-0723	740-0723	-	-
					.022" x .028"	741-0721	740-0721	-	-
 convertible	7 7	-14°	-	14°	.018" x .028"	741-0803	740-0803	741-0804	740-0804
					.022" x .028"	741-0801	740-0801	741-0802	740-0802

Mandibular

 convertible	6 6	-30°	-	4°	.018" x .028"	741-1703	740-1703	-	-
					.022" x .028"	741-1701	740-1701	-	-
 convertible, for Lip Bumper ø .045 inch / ø 1.15 mm	6 6	-30°	-	4°	.018" x .028"	741-1707	740-1707	-	-
					.022" x .028"	741-1705	740-1705	-	-
 convertible	7 7	-30°	-	4°	.018" x .028"	741-1803	740-1803	741-1804	740-1804
					.022" x .028"	741-1801	740-1801	741-1802	740-1802

Tulip buccal tubes



Maxillary	Tooth	Torque	Angulation	Distal offset	Slot	on large Pad		on Big Foot pad	
						right	left	v	left
 Tulip buccal tube	76 67	-14°	-	14°	.018" x .028"	743L0754	742L0754	743B0754	742B0754
					.022" x .028"	743L0752	742L0752	743B0752	742B0752