

ADHESIVES, CEMENTS & PRIMERS

Composites

ENLIGHT/ENLIGHT LV ORMCO (LIGHT-CURING ADHESIVES)


Ormco joined forces with Kerr R&D to develop Enlight and Enlight LV. These two materials contain a highly active initiator system for continuous curing in the absence of light for complete polymerization even with metal brackets.

Enlight

- ▶ Highly filled light-cure adhesive.
- ▶ Superb handling characteristics.
- ▶ Highly active initiator system.
- ▶ Unidose™ and syringe delivery options.

Enlight LV

- ▶ Flowable low-viscosity paste designed for lingual retainers and indirect bonding.
- ▶ Highly active initiator system.
- ▶ Unidose delivery for easy placement.



Designation	Kit Content	Part Number
Enlight syringe kit	4 syringes of paste (4 g) Ortho Solo (5 ml) Mixing pads and spatulas (75) Etching solution (10 ml) Applicator brushes (75)	740-0198
Enlight LV unidose kit Flowable paste (Low Viscosity)	6 boxes of adhesive paste (10 x 2 mg) Ortho Solo (5 ml) Mixing pads and spatulas (75) Etching solution (10 ml) Applicator brushes (75) Unidose dispensing gun	740-0299
Enlight unidose kit	6 boxes of adhesive paste (10 x 2 mg) Ortho Solo (5 ml) Mixing pads and spatulas (75) Etching solution (10 ml) Applicator brushes (75) Unidose dispensing gun	740-0199
Enlight introductory kit	1 syringe of adhesive paste (4 g) Ortho Solo (5 ml)	740-0179

Refill	Part Number
Enlight - Syringe (4g)	740-0195
Enlight - Unidose (20 x 2 mg)	740-0178
1 Ortho Solo bottle (5 ml)	740-0271
75 Mixing pads	740-0124
Bottle of etching solution (10 ml)	740-0038
Brushes (75)	740-0041
Unidose dispensing gun	740-0200
Enlight LV - Unidose (20 x 2 mg)	740-0278