








Memory Titanol® Spring for molar uprighting* (2nd Generation)

Flexible system for mesialization, intrusion or extrusion of molars.



Depending on the activation of the spring the appliance can achieve additional intrusive or extrusive moments.

| Order No. | | Content | Description |
|----------------------------------|---|---|---|
| Slot .018" 307-1011 |  |  | Set with cross tube |
| Slot .022" 307-1013 | | | |
| Slot .018" 307-1010 |  | 10 | Uprighting spring |
| Slot .022" 307-1012 | | | |
| Slot .018" 760-0062 |  | 10 | Cross tube |
| Slot .022" 760-0063 | | | |
| .016" x .016" 307-1009 |  | 10 | Molar uprighting spring (.016" x .016" – Stiff) Highly resilient, for the auxilliary slot of the Quick brackets |

*acc. to Prof. Dr. F. G. Sander