


Micro Sprint® Brackets



Info
N° 219

Roth Prescription

Farbmarkierung | colour-marks



Micro Sprint® Brackets Roth Prescription

OK / Maxillary					Slot .018"		Slot .022"		
Tooth		Torque	Angulation	In/Out	Rotation	Order No.		Order No.	
						Right	Left	Right	Left
1	Centrals	+12°	+5°	0,7	—	784-0103	783-0103	784-0101	783-0101
2	Laterals	+8°	+9°	1,0	—	784-0203	783-0203	784-0201	783-0201
3	Cuspids	-2°	+11°	0,5	4° mesial	784-0303	783-0303	784-0301	783-0301
	Cuspids + hook	-2°	+11°	0,5	4° mesial	784-0313	783-0313	784-0311	783-0311
4	Bicuspids	-7°	0°	0,5	2° distal	784-0503	783-0503	784-0501	783-0501
	Bicuspids + hook	-7°	0°	0,5	2° distal	784-0513	783-0513	784-0511	783-0511
5	Bicuspids	-7°	0°	0,5	2° distal	784-0503	783-0503	784-0501	783-0501
	Bicuspids + hook	-7°	0°	0,5	2° distal	784-0513	783-0513	784-0511	783-0511

UK / Mandibular					Slot .018"		Slot .022"		
Tooth		Torque	Angulation	In/Out	Rotation	Order No.		Order No.	
						Right	Left	Right	Left
1	Centrals	0°	0°	1,2	—	783-1203	783-1203	783-1201	783-1201
2	Laterals	0°	0°	1,2	—	783-1203	783-1203	783-1201	783-1201
3	Cuspids	-11°	+7°	0,7	2° mesial	784-1403	783-1403	784-1401	783-1401
	Cuspids + hook	-11°	+7°	0,7	2° mesial	784-1413	783-1413	784-1411	783-1411
4	Bicuspids	-17°	0°	0,7	4° distal	784-1503	783-1503	784-1501	783-1501
	Bicuspids + hook	-17°	0°	0,7	4° distal	784-1513	783-1513	784-1511	783-1511
5	Bicuspids	-22°	0°	0,7	4° distal	784-1603	783-1603	784-1601	783-1601
	Bicuspids + hook	-22°	0°	0,7	4° distal	784-1613	783-1613	784-1611	783-1611


Fälle | Cases

Variation	Slot .018"			Slot .022"		
Cases	1	5	10	1	5	10
1	706-1094	706-1095	706-1096	706-1103	706-1104	706-1105
3	706-1097	706-1098	706-1099	706-1106	706-1107	706-1108
3-5	706-1100	706-1101	706-1102	706-1109	706-1110	706-1111



McLaughlin/Bennett/Trevisi Prescription

Farbmarkierung | colour-marks



Micro Sprint® Brackets McLaughlin/Bennett/Trevisi Prescription

OK / Maxillary					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order No.		Order No.	
					Right	Left	Right	Left
1 Centrals	+17°	+4°	0,8	—	778-0103	777-0103	778-0101	777-0101
2 Laterals	+10°	+8°	1,0	—	778-0203	777-0203	778-0201	777-0201
3	Cuspids	-7°	+8°	0,5	778-0303	777-0303	778-0301	777-0301
	Cuspids + hook	-7°	+8°	0,5	778-0313	777-0313	778-0311	777-0311
4	Bicuspid	-7°	0°	0,5	778-0403	777-0403	778-0401	777-0401
	Bicuspid + hook	-7°	0°	0,5	778-0413	777-0413	778-0411	777-0411
5	Bicuspid	-7°	0°	0,9	778-0503	777-0503	778-0501	777-0501
	Bicuspid + hook	-7°	0°	0,9	778-0513	777-0513	778-0511	777-0511

UK / Mandibular					Slot .018"		Slot .022"	
Tooth	Torque	Angulation	In/Out	Rotation	Order No.		Order No.	
					Right	Left	Right	Left
1 Centrals	-6°	0°	1,2	—	777-1203	777-1203	777-1201	777-1201
2 Laterals	-6°	0°	1,2	—	777-1203	777-1203	777-1201	777-1201
3	Cuspids	-6°	+3°	0,7	778-1403	777-1403	778-1401	777-1401
	Cuspids + hook	-6°	+3°	0,7	778-1413	777-1413	778-1411	777-1411
4	Bicuspid	-12°	+2°	0,7	778-1503	777-1503	778-1501	777-1501
	Bicuspid + hook	-12°	+2°	0,7	778-1513	777-1513	778-1511	777-1511
5	Bicuspid	-17°	+2°	0,7	778-1603	777-1603	778-1601	777-1601
	Bicuspid + hook	-17°	+2°	0,7	778-1613	777-1613	778-1611	777-1611

Fälle | Cases

Variation	Slot .018"			Slot .022"		
Cases	1	5	10	1	5	10
	706-1112	706-1113	706-1114	706-1121	706-1122	706-1123
3	706-1115	706-1116	706-1117	706-1124	706-1125	706-1126
3-5	706-1118	706-1119	706-1120	706-1127	706-1128	706-1129