












Buccal tubes - special execution





Maxillary	Tooth	Torque	Angulation	Distal offset	Slot	for welding		on small pad	
						right	left	right	left
 convertible	6 6	-10°	-	8°	.018" x .028"	743-0703	742-0703	-	-
					.022" x .028"	743-0701	742-0701	-	-
 occlusal headgear tube convertible	6 6	0°	-	7°	.022" x .028"	759-1320	759-1220	-	-
	76 67	-10°	-	7°	.018" x .028"	-	-	757-8312	757-8212
 Mini buccal tubes	7 7	0°	-	0°	.018" x .028"	-	-	748-9311	748-9211
					.022" x .028"	-	-	748-9321	748-9221
 Comfort Line	7 7	-9°	-	10°	.022" x .028"	-	-	748-8322	748-8222

Mandibular

 convertible	6 6	-25°	-	6°	.018" x .028"	743-1703	742-1703	-	-
					.022" x .028"	743-1701	742-1701	-	-
 convertible	6 6	0°	-	7°	.018" x .028"	756-1310	756-1210	-	-
					.022" x .028"	756-1320	756-1220	-	-
 Mini buccal tubes	7 7	0°	-	0°	.018" x .028"	-	-	748-9311	748-9211
					.022" x .028"	-	-	748-9321	748-9221
 Comfort Line	7 7	-22°	-	12°	.018" x .028"	-	-	758-8314	758-8214

Tulip buccal tubes



Maxillary	Tooth	Torque	Angulation	Distal offset	Slot	on large pad		on Big Foot pad	
						right	left	right	left
	76 67	-10°	-	0°	.018" x .028"	743L0724	742L0724	743B0724	742B0724
					.022" x .028"	743L0722	742L0722	743B0722	742B0722
	76 67	-10°	-	6°	.018" x .028"	743L0734	742L0734	743B0734	742B0734
					.022" x .028"	743L0732	742L0732	743B0732	742B0732