

ADHESIVES, CEMENTS & PRIMERS

Primers

ORTHO SOLO AND ETCHING SOLUTIONS

Developed in partnership with the Kerr R&D department, Ortho Solo™ is a universal sealant and bond enhancer recommended for buccal and lingual bonding with chemical-cure systems by mixing or light-curing. It is used after etching.

- ▶ Fluoride release to prevent decalcification.
- ▶ Displaces moisture and does not form droplets.
- ▶ Perfect viscosity for easy application and to minimize bracket drift.
- ▶ Glass filler for optimal adhesion.
- ▶ No light-curing required with chemical adhesives.



Designation	Kit Content	Part Number
Complete kit	4 bottles (5 ml) Disposable wells (25) Applicator brushes (75)	740-0270
Ortho Solo sticks	Reduces the risk of cross-contamination Each Ortho Solo stick contains enough sealant/ bond enhancer for one arch 10 sticks (1,2 ml each)	740-0269
Ormco gel etch	Kit 2 syringes (2 x 5 ml / 25 tips)	740-0150
Ormco gel etch	Tips (25)	740-0134

Refill	Part Number
1 Ortho Solo bottle (5 ml)	740-0271