

TUBES

Titanium Tubes

TITANIUM TUBES

Titanium tubes are as resistant as stainless steel but twice as resilient. In addition, titanium is biocompatible avoiding allergies associated with nickel. Adhesion of bonding materials is enhanced due to the wettability of titanium.

- ▶ Wider base for 1st molar tubes.
- ▶ Teardrop-shaped mandibular tube to reduce occlusal interference.
- ▶ Funneled opening for easy archwire insertion.
- ▶ Notch for easy handling of the tube.
- ▶ 1st and 2nd molars available in .018" x .025" and .022" x .028".

1st upper molar



1st lower molar



2nd upper molar



2nd lower molar

