

Titanol® Low Force

Einzel verpackt / individually packed




Beschreibung description	Straight-Arch- Form		Euro-Smile- Form		Profil profile	ø mm	ø inch	Stück content
	Order No. OK / Maxillary	Order No. UK / Mandibular	Order No. OK / Maxillary	Order No. UK / Mandibular				
Titanol® Low Force mit 30 % weniger Rückstellkraft als Titanol Superelastic bis zu 30% weniger Reibung Titanol® Low Force with 30 % less restoring force than Titanol Superelastic up to 30 % less friction	280-0840	280-0940	-	-		0,40	.016"	10
	280-1440	280-1540	281-1440	281-1540		0,41 x 0,41	.016" x .016"	10
	280-2040	280-2140	281-2040	281-2140		0,41 x 0,56	.016" x .022"	10
	280-2046	280-2146	-	-		0,46 x 0,64	.018" x .025"	10

Titanol® Triple Force

Straight-Arch-Form


Einzel verpackt / individually packed



Die 3 Kraftzonen / the 3 force zones

Mit den drei Kraftsegmenten für eine biologisch günstige Kraft auf die Zähne. Light für die Front, medium für den Seitenzahnbereich und heavy für den posterioren Teil.

With the three force segments for a biologically favorable force on the teeth. Light for anteriors, medium for bicuspid and heavy for the posterior segment.



Beschreibung description	Order No. OK / Maxillary	Order No. UK / Mandibular	Profil profile	ø mm	ø inch	Stück content
Triple Force Drähte in Straight-Arch-Bogenform	276-1440	276-1540		0,41 x 0,41	.016" x .016"	10
Triple Force arch wires in Straight-Arch form	276-2041	276-2141		0,41 x 0,56	.016" x .022"	10

* bis zu 30% weniger Reibung / up to 30 % less friction