

# Acrylics technique

Learn more:



## STEADY-RESIN S

Cold-curing resin for the "salt and pepper" technique. This polymer meets the requirements of the "salt and pepper" technique, such as high stability and optimum wettability.

// polymer

8134.1  
1 kg

CE



## STEADY-RESIN M

Cold-curing resin for the "premix" technique. Characterised by a wide processing range and excellent modelling properties.

// polymer

8135.1  
1 kg

CE



## STEADY-RESIN VARIO

Cold-curing resin for orthodontics with prolonged processing time, can be used both for "salt and pepper" and for the modelling techniques.

// polymer

8140.1  
1 kg

CE



## STEADY-RESIN S+M

Universal monomer, can be processed with all STEADY-RESIN polymer versions. The clear monomer can be mixed with STEADY-RESIN coloured monomers.

// monomer, clear

8136.1  
500 ml

// monomer, pink

8137.1  
500 ml

CE