


## BioStarter®, BioTorque®


Einzeln verpackt / individually packed



Beschreibung description	Order No. OK / Maxillary	Order No. UK / Mandibular	Order No. OK / Maxillary	Order No. UK / Mandibular	Profil profile	Kraft force	ø inch	Stück Content
<b>BioStarter®</b>	203-0825	203-0925	203-1825	203-1925	●	20 g	.010"	10
	203-0830	203-0930	203-1830	203-1930	●	30 g	.012"	10
	203-0835	203-0935	203-1835	203-1935	●	40 g	.014"	10
	203-0840	203-0940	203-1840	203-1940	●	50 g	.016"	10
<b>BioTorque®</b>	203-2044	203-2144	203-4044	203-4144	■	140 g	.017"x.025"	10
	203-2048	203-2148	203-4048	203-4148	■	170 g	.019"x.025"	10

## Stahlbögen

Code:



Beschreibung description	Order No. OK / Maxillary	Order No. UK / Mandibular	Order No. OK / Maxillary	Order No. UK / Mandibular	Profil profile	ø mm	ø inch	Stück content
<b>Stahlbögen Steel arch wires</b>	202-3748	202-3848	201-3748	201-3848	■	0,48 x 0,64	.019"x.025"	10