


# Tulip buccal tubes

Info  
N° 208

## Tulip buccal tubes

The Tulip buccal tubes owe their name to the tulip-shaped entrance funnel. It has been dimensioned especially large to allow quick and easy insertion of the archwires. Your patients will be delighted by the extremely flat design which guarantees a high level of comfort.



Maxillary	Tooth	Torque	Angulation	Distal offset	Slot	on large pad		on Big Foot pad	
						right	left	right	left
						76   67	0°	-	0°
				.022" x .028"	743L0712	742L0712	743B0712	742B0712	
76   67	-10°	-	0°	.018" x .028"	743L0724	742L0724	743B0724	742B0724	
				.022" x .028"	743L0722	742L0722	743B0722	742B0722	
76   67	-10°	-	6°	.018" x .028"	743L0734	742L0734	743B0734	742B0734	
				.022" x .028"	743L0732	742L0732	743B0732	742B0732	
76   67	-14°	-	10°	.018" x .028"	743L0744	742L0744	743B0744	742B0744	
				.022" x .028"	743L0742	742L0742	743B0742	742B0742	
76   67	-14°	-	14°	.018" x .028"	743L0754	742L0754	743B0754	742B0754	
				.022" x .028"	743L0752	742L0752	743B0752	742B0752	

## Mandibular

7   7	-10°	-	0°	.018" x .028"	742T0724	743T0724	742B0724	743B0724
				.022" x .028"	742T0722	743T0722	742B0722	743B0722
76   67	-20°	-	0°	.018" x .028"	743L1714	742L1714	743B1714	742B1714
				.022" x .028"	743L1712	742L1712	743B1712	742B1712
76   67	-25°	-	6°	.018" x .028"	743L1724	742L1724	743B1724	742B1724
				.022" x .028"	743L1722	742L1722	743B1722	742B1722
76   67	-30°	-	6°	.018" x .028"	743L1734	742L1734	743B1734	742B1734
				.022" x .028"	743L1732	742L1732	743B1732	742B1732