


BANDS

Ormco Bands

WASHBON™ FIRST MOLARS

STANDARD SELECTION
Each kit contains an equal number of left bands and right bands.

- ▶ Shaped to tooth anatomy for a snug, secure fit and optimal retention.
- ▶ Right and left anatomy is proportioned in all sizes for optimum adjustment.
- ▶ Perfect reproduction of the three buccal cusps of the 1st lower molar on the band.



1st maxillary and mandibular molars (suggested kits)

Kit/Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
100 bands	2	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	6	6	6	6	4	4	4	4	2	2	2	2	2	2	2
200 bands	2	2	2	2	2	2	4	4	6	6	6	8	10	10	10	10	12	12	12	10	10	10	8	8	6	6	6	4	4	2	2	2
400 bands	2	2	4	4	6	6	8	8	10	14	14	16	20	20	24	24	26	26	24	24	20	20	18	14	12	10	8	6	4	2	2	2
800 bands	4	4	6	6	10	10	12	18	24	28	32	34	40	44	48	50	50	50	50	50	44	40	36	26	24	20	14	10	4	4	4	4
100 bands												6	8	8	8	10	10	10	10	8	8	8	6									

Size	Side	Tooth	U/L	Part Number*
01 to 32	Left	1st molar	U6	378-00XX
01 to 32	Right	1st molar	U6	377-00XX
01 to 32	Left	1st molar	L6	374-00XX
01 to 32	Right	1st molar	L6	373-00XX

* The two "XX" at the end of the part number represent the size of the band.


BANDS

Ormco Bands

WASHBON™ SECOND MOLARS

STANDARD SELECTION
Each kit contains an equal number of left bands and right bands.

► Occlusal festoon at distal portion of upper 2nd molar band helps eliminate metal overhang and occlusal interference.



2nd maxillary and mandibular molars (suggested kits)

Kit/Number maxillary	97	98	99	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28
Kit/Number mandibular	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
100 bands	2	2	2	2	2	2	2	2	2	2	4	4	4	4	4	4	6	6	6	6	4	4	4	4	2	2	2	2	2	2	2	2
200 bands	2	2	2	2	2	2	4	4	6	6	6	8	10	10	10	10	12	12	12	10	10	10	8	8	6	6	6	4	4	2	2	2
400 bands	2	2	4	4	6	6	8	8	10	14	14	16	20	20	24	24	26	26	24	24	20	20	18	14	12	10	8	6	4	2	2	2
800 bands	4	4	6	6	10	10	12	18	24	28	32	34	40	44	48	50	50	50	50	50	44	40	36	26	24	20	14	10	4	4	4	4
100 bands											6	8	8	8	10	10	10	10	8	8	8	6										

Size	Side	Tooth	U/L	Part Number*
97 to 28	Left	2nd molar	U7	372-00XX
97 to 28	Right	2nd molar	U7	371-00XX
01 to 32	Left	2nd molar	L7	376-00XX
01 to 32	Right	2nd molar	L7	375-00XX

* The two "XX" at the end of the part number represent the size of the band.