

ZOO PACK ELASTICS





































Shaded sizes are available in Ormco Color Packs

Chipmunk	636-0113	Bear	636-0122
Rabbit	636-0114	Impala	636-0125
Fox	636-0115	Ram	636-0126
Penguin	636-0116	Moose	636-0127
Monkey	636-0117	Buffalo	636-0128

- High-quality surgical latex
- Consistent force pull
- Force ratings based on elastic being extended three times listed diameter
- Available in 10 sizes of assorted Color Packs (pink, orange, green and purple)
- Zoo Packs contain 50 bags of 100

INTRAORAL ELASTICS

EXTRAORAL ELASTICS

strength	light	medium	med / heavy	heavy	extra heavy	light	heavy
SIZE	2oz — 60g	3oz — 85g	3.5oz — 100g	4.5oz — 130g	6 oz — 170 g	8 oz — 230 g	14 oz — 400 g
1/8" 3.18 mm	 Hummingbird 630-0010		 Chipmunk 630-0030				
3/16" 4.76 mm	 Quail 630-0011	 Otter 630-0020	 Rabbit 630-0031	 Kangaroo 630-0040	 Impala 630-0050	 Cougar 635-0058	
1/4" 6.35 mm	 Owl 630-0012	 Seal 630-0021	 Fox 630-0032	 Bear 630-0041	 Ram 630-0051	 Leopard 635-0059	
5/16" 7.94 mm	 Parrot 630-0013	 Dolphin 630-0022	 Penguin 630-0033	 Zebra 630-0042	 Moose 630-0052	 Panther 635-0060	 Walrus 635-0065
3/8" 9.35 mm	 Road Runner 630-0014	 Turtle 630-0023	 Monkey 630-0034	 Camel 630-0043	 Buffalo 630-0053	 Tiger 635-0061	 Elephant 630-0066
1/2" 12.7 mm	 Peacock 630-0015		 Donkey 630-0035			 Lion 635-0062	 Whale 630-0067
5/8" 15.9 mm	 Eagle 630-0016		 Llama 630-0036				
3/4" 19.1 mm	 Ostrich 630-0017		 Giraffe 630-0037				

Animal Zoo Packs	Master Pouches of 50 Zoo Packs	Color Packs of 40 Zoo Packs
LIGHT - 2 OZ		
Hummingbird	630-0010	
Quail	630-0011	
Owl	630-0012	
Parrot	630-0013	
Road Runner	630-0014	
Peacock	630-0015	
Eagle	630-0016	
Ostrich	630-0017	
MEDIUM - 3 OZ		
Otter	630-0020	
Seal	630-0021	
Dolphin	630-0022	
Turtle	630-0023	
MEDIUM HEAVY - 3.5 OZ		
Chipmunk	630-0030	636-0013
Rabbit	630-0031	636-0014
Fox	630-0032	636-0015
Penguin	630-0033	636-0016
Monkey	630-0034	636-0017
Donkey	630-0035	
Llama	630-0036	
Giraffe	630-0037	

Animal Zoo Packs	Pouches of 15 Zoo Packs	Master Pouches of 50 Zoo Packs	Color Packs of 40 Zoo Packs
HEAVY - 4.5 OZ			
Kangaroo		630-0040	
Bear		630-0041	636-0122
Zebra		630-0042	
Camel		630-0043	
EXTRA HEAVY - 6 OZ			
Impala		630-0050	636-0125
Ram		630-0051	636-0126
Moose		630-0052	636-0127
Buffalo		630-0053	636-0128
LIGHT - 8 OZ HEADGEAR			
Cougar	635-0058	634-0000	
Leopard	635-0059	634-0001	
Panther	635-0060	634-0002	
Tiger	635-0061	634-0003	
Lion	635-0062	634-0004	
HEAVY - 14 OZ HEADGEAR			
Walrus	635-0065	634-0012	
Elephant	635-0066		
Whale	635-0067	634-0014	



ZOO PACK DISPENSER

- Convenient, clear plexi rack
- Holds 6 Zoo Pack boxes
- Easily mounts on wall or sits on counter top (Zoo Packs not included)

Dispenser (Zoo Packs not included)	630-0000
------------------------------------	----------

1
2
3
4
5
6
7
8
9
10
11
12
13